

TENSOMAX


We Measure Texture

Universal Texture Analyzer

Model: TexMax TX-100

Features:

- # Load capacity: 100kgf/1KN
- # Speed range: 0.1 — 1000 mm/min
- # Servo motor precise control
- # Force resolution up to 1/100,000
- # Professional yet easy software
- # Easy to use Windows-based software.
- # Ready to use international method library



Cat No. : 0005289






Software Capabilities:

Our unique TexMax software can measure the following items:

- * Hardness
- * Firmness
- * Springiness
- * Cohesiveness
- * Resilience
- * Fracturability
- * Adhesiveness
- * Gumminess
- * Chewiness
- * and much more properties.



Features:

- * Everybody, even with low knowledge of test theory, can do the tests.
- * Calibration of the force, extension and speed in the software is very convenient.
- * No need to engage for setup the methods. Since our engineers have made a lot of ready-to-use methods.
- * During every test, user can watch the test presentation via multi-graph.
- * A video camera can capture test videos during test, and software will catch them in the record section for future evaluation and history.
- * All test data and records will be saved and they can be sorted and filtered by the time, method and test name filters.
- * And the most important thing is: our software is professional yet easy and user friendly.


TENSOMAX




Cat No. : 0502280
Cat name: Back Extrusion



Cat No. : 0502281
Cat name: Spread-cone-45




Cat No. : 0502282
Cat name: cheese cutting



Cat No. : 0502283
Cat name: cheese compression



Cat No. : 0502284
Cat name: cheese cone spread



Cat No. : 0502286
Cat name: cheese penetration



Cat No. : 0502289
Cat name: Cheese Extensibility

TENSOMAX



Cat No. : 0502250
Cat name: Spread-cone-45



Cat No. : 0502251
Cat name: freshness



Cat No. : 0502252
Cat name: firmness and
springiness of cake




Cat No. : 0502253
Cat name: Biscuit Snap Strength



Cat No. : 0502253
Cat name: Chocolate Snap
Strength




Cat No. : 0502260
Cat name: Snack Crispness



Cat No. : 0502263
Cat name: Bloom hardness
Strength of bread to AACC




Cat No. : 0502264
Cat name: Bread V-squeeze




Cat No. : 0502268
Cat name: dough extensograph


TENSOMAX




Cat No. : 0502270
Cat name: Warner-Bratzler



Cat No. : 0502271
Cat name: Warner-Bratzler
_V slot



Cat No. : 0502274
Cat name: meat fracture.



Cat No. : 0502276
Cat name: Volodkevich Bite Jaws



Cat No. : 0502278
Cat name: egg compression

TENSO MAX




Cat No. : 0502290
Cat name: burger consistency



Cat No. : 0502283
Cat name: burger compression



Cat No. : 0502295
Cat name: burger quality
in teeth




Cat No. : 0502283
Cat name: Sandwich
Compression



Cat No. : 0502297
Cat name: pizza tensile

TENSOmax



Cat No. : 0502320
Cat name: kramer shear
mini cell



Cat No. : 0502321
Cat name: bending vertical




Cat No. : 0502322
Cat name: bending horizontal




Cat No. : 0502323
Cat name: Spaghetti tensile

TENSOMAX



Cat No. : 0502310
Cat name: flesh firmness



Cat No. : 0502311
Cat name: fruit hardness

TENSOmax



Cat No. : 0502341
Cat name: Kramer Shear Cell



Cat No. : 0502100
Cat name: Video recording module